

CureMed Nordic AB is a Swedish medtech company with head office and warehouse in Gothenburg. In our range, you will find world-leading products within healthcare and medical technology, many of which can be customized in production. Our goal is to always offer innovative products from leading manufacturers in a cost-effective manner to the healthcare sector.

# Replacement Clamp™

- Replacement clamp is a repair clamp for dialysis catheters\*.
- Replacement clamp makes it possible to replace a broken clamp on the dialysis catheter without removing the luer lock connector.
- New ergonomic design.
- Manufactured in Sweden.
- CE certified.



**Step 1:**

Remove the old clamp. Slide in the Replacement clamp from the side.

**Step 2:**

Fix the tube in the notch of the Replacement clamp. The clamp is now ready to use.



*\*Replacement clamp is for temporary use only. A damaged catheter should always be replaced. Follow the instructions of your catheter supplier carefully.*

**ORDERING INFORMATION**

Product code	Description	Packaging
RCR01	Replacement clamp, red	2 pcs / bag
RCB02	Replacement clamp, blue	2 pcs / bag
RCM03	Replacement clamp, red & blue	2 pcs (1 red & 1 blue) / bag