

CureMed Nordic AB is a Swedish medtech company with head office and warehouse in Gothenburg. In our range, you will find world-leading products within healthcare and medical technology, many of which can be customized in production. Our goal is to always offer innovative products from leading manufacturers in a cost-effective manner to the healthcare sector.

Replacement Clamp™

- Replacement clamp is a repair clamp for dialysis catheters*.
- Replacement clamp makes it possible to replace a broken clamp on the dialysis catheter without removing the luer lock connector.
- New ergonomic design.
- Manufactured in Sweden.


CE FDA reg number
3015532181

Step 1:

Remove the old clamp.
Slide in the Replacement clamp from the side.

Step 2:

Fix the tube in the notch of the Replacement clamp.
The clamp is now ready to use.



**Replacement clamp is for temporary use only. A damaged catheter should always be replaced. Follow the instructions of your catheter supplier carefully.*

ORDERING INFORMATION

Product code	Description	Packaging
RCM03	Replacement clamp, red & blue	2 pcs (1 red & 1 blue) / bag